

Doric Gingival Retraction Paste

Instructions for Use

PRODUCT DESCRIPTION

Doric Gingival Retraction Paste is a retraction agent for sulcus widening. Doric Gingival Retraction Paste has an astringent effect on the gingiva.

INDICATIONS

Dental treatment cases that require final impressions of prepared teeth e.g. after tooth loss.

CONTRA-INDICATIONS

Do not use the material where there are contact allergies or known allergies to any of the components.

INTENDED PURPOSE

For temporary marginal gingival retraction and provision of a dry sulcus environment in the case of healthy periodontium.

LIMITATIONS OF USE

Do not use in the case of exposed bone, advanced periodontopathy or exposed furcation defects.

INTENDED USERS

Dentist.

PATIENT TARGET GROUP

Persons receiving treatment as part of a dental procedure.

RECOMMENDED USE

1. Prepare the tooth as required.
2. Place the compule in a suitable dispenser.
3. Bend the tip of the compule into the required position if necessary.
4. Dispense a small amount of paste onto a mixing block and discard.
5. Optional: If needed for the clinical situation, place a retraction cord in the sulcus.
6. Rinse the sulcus with water and then dry the tissue.
7. Place the paste slowly around the preparation in the sulcus until the sulcus is completely filled. Where necessary,

use compression caps or other compression products.

8. Leave the paste in place for approx. 02:00 minutes, to enable the gingival retraction. While it is setting, the area being treated must be kept dry.

Caution: Leaving the paste in the mouth for too long can lead to irritations of the gingiva. Do not exceed the intraoral exposure time.

9. Remove the paste completely from the sulcus and mouth area with a mixture of air and water and by using suction.

Note: Completely remove the paste.

10. Residues of the paste can lead to inaccuracies or delays in the setting of the impression material. Dry the gingiva with air and check the retraction result. If a stronger retraction is required, repeat the process from step 7. If necessary, a thin retraction cord can be used in addition.
11. Taking the impression.

TIMING

Reaction time: 02:00 minutes

NOTES FOR USE

Use the compule immediately after opening the aluminum bag. Discard any residual material in the compule.

INTERACTIONS

Residues of the paste can lead to delays in the setting of the impression material. Remove the paste completely from the sulcus and mouth area with a mixture of air and water and by using suction.

RESIDUAL RISKS / SIDE EFFECTS

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

WARNINGS / PRECAUTIONS

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

PERFORMANCE CHARACTERISTICS OF THE DEVICE

Serves to widen the sulcus and has an astringent effect on the gingiva.

STORAGE / DISPOSAL

- Store in a dry place at room temperature (15 - 25 °C/59 - 77 °F)! Do not use after the expiration date.
- Use the compule immediately after opening the aluminum bag. Discard any residual material in the compule.
- Disposal must comply with national regulations.

COMPOSITION

Kaolin clay, water, aluminum chloride hexahydrate, additives.

PACKAGING

REF: 0546 20 compules @ 0.37 g

Doric is an internationally registered trademark of Davis Schottlander & Davis Ltd.

