Davis Schottlander & Davis Ltd.
Letchworth, Herts, UK
Tel: +44 (0)1462 480848
www.schottlander.com

schottlander

Pegasus Plus Repair Acrylic Instructions for Use

INTRODUCTION

Pegasus Plus repair acrylic from Schottlander is a repair acrylic based on methyl methacrylate and complies with EN ISO 20795-1:2013, Type 2 : Class 1. This product is specifically formulated for use in dentistry and for use by dental care professionals only.

INTENDED USE

To repair part of a denture (partial or complete) that is used to replace missing teeth in humans.

INDICATIONS FOR USE

Adults or children requiring the introduction of teeth in the oral cavity for functional reasons and where a denture base is required.

CONTRA-INDICATIONS

People with known allergies to methacrylates should only use with suitable protection. Certain patients are unable to tolerate the wearing of appliances made from methacrylates.

DANGER

This product contains methyl methacrylate and the liquid is **highly flammable** and classed as **irritant**. Read health & safety section of these directions and request a copy of the Safety Data Sheet before use.

PREPARATION

Assemble the broken parts and cast the model in type 3 plaster. Cut back the broken area to leave a 1.5mm gap and bevel the edges to form a shallow V. Protect the fitting surfaces and nearby teeth with a thin film of petroleum jelly. Coat the model with Schottlander Isolating Solution and secure the denture to the model with sticky wax.

MIXING RATIO

13g powder: 10ml liquid.

MIXING METHOD A

Moisten the contact surfaces of the repair with Pegasus Plus repair liquid. Add Pegasus Plus repair powder and continue to moisten and add the powder until the repair has been built up in excess of the thickness needed.

MIXING METHOD B

Pour 10 drops of liquid into a dry Dappen glass. Add powder evenly until all monomer is absorbed and there is a dry excess. Reverse the Dappen glass and tap in the palm of the hand to remove excess powder. Add four drops of monomer and stir gently for 20-30 seconds with spatula. Avoid trapping air. The approximate working time is 1 minute 30 seconds to 2 minutes.

APPLICATION

The mix should be free flowing - if sluggish, discard and make a new mix. Apply using a thin bladed spatula and slightly overbuild. Once in position do not disturb. Excess of flow can be controlled with a dusting of dry powder.

CURING IN FLASK

Best results are obtained with flask curing. Immerse into water at 40°C bring the pressure to 50psi (3.4 bar or 3.5kg/cm²). Hold for about 10 minutes. This may vary according to the type of flask used.

AIR CURING

It is also possible to obtain a porosity free repair by leaving undisturbed to cure in the air. In common with other acrylics the setting time is sensitive to temperature and at 23°C is approximately 15 minutes. Warm conditions decrease this time and vice versa.

FINISHING

When completely set and hard remove from model and trim and polish in the usual way.

STORAGE

Store in a dry place at room temperature (max 30°C) and avoid prolonged exposure to sunlight. Keep containers closed when not in use.

HEALTH & SAFETY

Warnings: Liquid contains methyl methacrylate monomer. H225 Flammable liquids. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation.

Precautions: P210: Keep away from heat/sparks/open flames/hot surfaces. No smoking. P261: Avoid breathing dust/fume/gas/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P303+361+353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P501: Dispose of contents/container in accordance with local regulation.

LOT NUMBERS

The lot number and the expiry date are shown on all containers.

PRODUCT CODES

POWDER: 500g

REF 0521-P Pink

REF 0521-V Pink Veined

REF 0521-TV Translucent Veined

REF 0521-C Clear

POWDER: 1kg

REF 0522-P Pink

REF 0522-V Pink Veined

REF 0522-TV Translucent Veined

REF 0522-C Clear

LIQUID

REF 0526 200ml

REF 0527 1 litre

REF 0529 2.5 litres

Pegasus is an internationally registered trade mark of Davis Schottlander & Davis Limited.

Where this product is shown as having been certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, this is exclusively for the indication(s) shown in the above Instructions for Use. Other non-medical uses ascribed to this device are not within the scope of CE certification, and users should be aware that product performance and/or safety has not been evaluated by SGS for those purposes.

Davis Schottlander & Davis Ltd.
Letchworth, Herts, UK
Tel: +44 (0)1462 480848
www.schottlander.com

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Pegasus är ett internationellt registrerat varumärke som tillhör Davis Schottlander & Davis Limited.
EC REP Schottlander (IRL) Limited. Suite 11, Parklands Office Park, Southern Cross Road, Bray, Co Wicklow, Eire.