Davis Schottlander & Davis Ltd.
Letchworth, Herts, UK
Tel: +44 (0)1462 480848
www.schottlander.com

schottlander

Pegasus Plus Denture Base Material Instructions for Use

INTRODUCTION

Pegasus Plus denture base material from Schottlander is a fast curing heat cured acrylic based on methyl methacrylate and complies with EN ISO 20795-1:2013, Type 1 : Class 1. This product is specifically formulated for use in dentistry and for use by dental care professionals only.

INTENDED USE

To form part of a denture (partial or complete) that is used to replace missing teeth in humans.

INDICATIONS FOR USE

Adults or children requiring the introduction of teeth in the oral cavity for functional reasons and where a denture base is required.

CONTRA-INDICATIONS

People with known allergies to methacrylates should only use with suitable protection. Certain patients are unable to tolerate the wearing of appliances made from methacrylates.

DANGER

This product contains methyl methacrylate and the liquid is **highly flammable** and classed as irritant. Read health & safety section of these instructions and request a copy of the Safety Data Sheet before use.

PHYSICAL DATA & PROCESSING TIMES

Powder:liquid ratio	2.34g per ml
Dough/Initial Packing time at 23°C	6 ± 2 mins
Working/Final Packing time at 23±2°C	20 mins

WAXING & INVESTING

Invest in gypsum and follow procedure for gypsum models. Use a proprietary denture flask.

WAXING REMOVAL

Immerse the flask in boiling water for five minutes. Remove and open the flask, lift out wax and flush thoroughly with boiling water to which detergent solution has been added, followed by clean boiling water. Allow to cool.

APPLY SCHOTTLANDER ISOLATING SOLUTION

Use a clean brush to apply one coat of Schottlander Isolating Solution to all plaster surfaces. Absorb excess from around teeth while it remains fluid. Make quite sure there is none on the exposed surfaces of the teeth. Allow to dry before packing.

MIXING

Caution: Liquid flash point 9.5°C - do not mix or store near naked flame. Recommended ambient temperature for mixing and packing is 23°C. Pour liquid into a mixing vessel and add appropriate powder in a slow stream (taking 30-40 seconds) until an excess appears on the surface. Hold the vessel in the hand, tap 3 or 4 times to bring excess liquid to the surface and add sufficient powder to absorb this liquid. Invert the vessel and discard any surplus powder. Thoroughly spatulate the mix for 1 minute. Cover the vessel and wait for the mix to reach the correct stage for packing and pressing (dough time). This stage is reached when the mix can be separated cleanly from the walls of the mixing vessel with a spatula.

PACKING

Care: Do not pack before the dough stage has been reached. The flask temperature should be 23°C.

Pack the material into the mould and cover with a thin polythene sheet. Assemble the flask and close slowly under a press. Open the flask, paint the model half of the flask with Schottlander Isolating Solution, remove separating sheet and cut away the flash. Then reassemble the flask, close under press and transfer to a security flask clamp. As an extra safeguard against porosity place a small additional amount of dough centrally in the flask. Do not pack after the end of the working time.

EXTENDED WORKING TIME

If storing Pegasus Plus denture base for extended use, ensure that dough stage is reached prior to refrigeration.

RECOMMENDED CURES

Normal dentures: Immerse the clamped flask in cold, warm or boiling water. Bring to the boil in not less than 5 minutes and simmer for a further 25 minutes. Residual monomer content < 0.6%.

Thick dentures: For thicker dentures place the clamped flask in boiling water and remove heat source for 20 minutes. Bring to the boil in not less than 5 minutes and continue boiling for a further 10 minutes. Residual monomer content < 0.6%.

Dentures with reduced residual monomer: Place the clamped flask in cold water and raise to 70°C. Hold this temperature for 90 minutes and then bring to the boil and continue boiling for 2 hours. Residual monomer content < 0.5%.

ALTERNATIVE CURES

Other longer cure techniques can, if desired, be employed with Pegasus Plus denture base material. When dry curing follow the instructions of the bath manufacturer. **Note:** After all curing cycles bench cool the flask for 30 minutes prior to immersing in cold water before deflasking. Finish and polish in the usual manner. Store finished denture wrapped in damp tissue in a polythene bag.

INSTRUCTIONS FOR ADDITION & REPAIRS

To obtain maximum repair strength, it is recommended that repairs or additions be made using Pegasus Plus heat cure denture base. An autopolymerising resin such as Pegasus Plus Repair Acrylic may also be used for this purpose.

STORAGE

Store at room temperature (max 30°C) and avoid prolonged exposure to sunlight. Keep containers tightly closed when not in use.

HEALTH & SAFETY

Hazards: Liquid contains methyl methacrylate monomer. H225 Highly flammable liquid or vapour.

H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation.

Precautions: P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking. P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P281 Use personal protective equipment as required. P303 + 361 + 353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P312 Call a POISON CENTRE or doctor/physician if you feel unwell. P330 Rinse mouth. P501 – Dispose of contents/container in accordance with local regulation.

LOT NUMBERS

The lot number and the expiry date are shown on all containers.

Pegasus is an internationally registered trade mark of Davis Schottlander & Davis Limited.

Where this product is shown as having been certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, this is exclusively for the indication(s) shown in the above Instructions for Use. Other non-medical uses ascribed to this device are not within the scope of CE certification, and users should be aware that product performance and/or safety has not been evaluated by SGS for those purposes.

EC REP Schottlander (IRL) Limited. Suite 11, Parklands Office Park, Southern Cross Road, Bray, Co Wicklow, Eire.

Davis Schottlander & Davis Ltd.
Letchworth, Herts, UK
Tel: +44 (0)1462 480848
www.schottlander.com

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Pegasus är ett internationellt registrerat varumärke som tillhör Davis Schottlander & Davis Limited.