enliven MC

Instructions for Use

Indications: enliven MC is a veneering ceramic porcelain for bonding to metal crowns and bridges.

Contra Indications: Do not use without a metal sub-frame. Not suitable for use with low fusing porcelains, PJC materials or with some non-precious alloys. Do not use with a sub-frame other than described in 1. Alloys below.

1. **Alloys:** Non Precious, Semi Precious & Precious Alloys with coefficient of expansion between $13.8 - 14.9 \times 10^{-6} \text{ }^\circ\text{C}^{-1}$ may be used. Prepare your copings and where indicated degas to the alloy manufacturer's directions for use.

2. **Opaque Applications:** Stir contents in the tub before each application. Apply an even thin layer completely masking the coping. When fired the 1st Opaque is shiny with a partial opaque coverage to the applied surface.


2nd Opaque: Apply a thin even layer masking the 1st Opaque. When fired the Opaque has a slight shine.

3. Dentine & Incisal Build Up

Build the dentine to full contour and reduce for incisal layer:

Figure 1 - when dentine thickness 0.5 to 0.7mm

Figure 2 - when dentine thickness >0.7mm

4. Small Corrections & Glaze may be fired at the same time.

Furnace temperatures vary. Adjust high temperature to give appearance in firing chart.

Oxidation	Refer to the manufacture's instructions for the alloy						
	Start Temp	Minimum Drying Time	Temp Rise $^\circ\text{C}/\text{min}$	Vacuum $^\circ\text{C}/\text{min}$ (720mm Hg)	High Temp	Hold Time Air	Appearance
1st Opaque Firing	450 $^\circ\text{C}$	8 minutes	80 $^\circ\text{C}$	Yes	980 $^\circ\text{C}$	1 minute	Shiny
2nd Opaque Firing	450 $^\circ\text{C}$	6 minutes	60 $^\circ\text{C}$	Yes	950 $^\circ\text{C}$	1 minute	Slight Shine
All other firings*	450 $^\circ\text{C}$	6 minutes	60 $^\circ\text{C}$	Yes	930 $^\circ\text{C}$	1 minute	Slight Shine

* Through Unifire® Technology. Precious alloys generally require the opaque temperatures to be fired 20 $^\circ\text{C}$ lower.

5. Shoulder Application

Isolate the gypsum die with the Schottlander ceramic separator pen using two coats. Once soaked in and dry, apply a very fine bead of enliven Shoulder Opaque Margin and fire at 930 $^\circ\text{C}$.

If a shoulder is required mix one part enliven MC Shoulder Stabiliser with four parts of the dentine shade required. Build up the dentine and enamel over the top and fire in one firing.

Note. If a higher proportion of Shoulder Porcelain is used then fire separately at 940-950 $^\circ\text{C}$.

6. **enliven Universal stains** (suitable for both metal and all ceramic systems) enliven stains are designed so that they may be used as powder or paste stains according to choice. For powder stains mix with enliven Glaze & Stain liquid. For paste stains mix with enliven Paste Stain liquid.

Coefficient of expansion between 25 $^\circ\text{C}$ – 500 $^\circ\text{C}$

enliven paste opaques, dentines, enamels, translucents, opal halo $13.4 \times 10^{-6} \text{ K}^{-1}$
Shoulder porcelains $13.8 \times 10^{-6} \text{ K}^{-1}$

Glass transition temperature

enliven paste opaques and powders 590 $^\circ\text{C}$

STORAGE

Store enliven powders between 4 $^\circ\text{C}$ and 29 $^\circ\text{C}$ and paste opaques between 15 $^\circ\text{C}$ and 23 $^\circ\text{C}$. Immediately after each usage replace all caps.

HEALTH AND SAFETY

Keep away from eyes since the granular nature of the powder or paste may cause irritation.

LOT NUMBERS

The lot number is shown on all packs.

This product is specifically formulated for use in dentistry.

enliven is an internationally registered trade mark of Davis Schottlander & Davis Limited.

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