

Doric Pre-cut impression separation wafers

INTRODUCTION

Doric Pre-cut impression separation wafers are made from 1mm polyethylene foam and pre-cut to size so that trimming is only carried out when required. They are manufactured to create an even space in the initial Doric putty impression. This guarantees an equal distribution of the Doric wash impression material and provides a high level of sub gingival detail.

Doric Pre-cut impression separation wafers are suitable for both upper and lower arch impressions. For the lower impression a pre-cut flap is easily removed.

CONTRA INDICATIONS

Patients with known allergy to polyethylene.
See Health & Safety below.

DIRECTIONS FOR USE

- Prior to preparation of the teeth select a suitable sized rigid impression tray. Check the size of the wafer against the tray and, if necessary, trim to shape. Allow for the rear 2mm of the tray to remain uncovered and ensure a complete impression of the teeth at the end of the arch. This will form a stop for any excess Doric wash material when taking the final impression. For lowers remove the central flap.
- Mix the Doric putty and distribute uniformly over the tray. Place a separation wafer in position and lightly press down so that it sticks to the surface. Insert the tray into the mouth and take the first impression. Remove the impression immediately the Doric putty starts to set.
- Remove the separation wafer from the impression. The surface will appear textured and this will offer increased surface area for retention of the second impression material. Check the fit of the tray in the mouth and remove any excess material with a scalpel.
- Clean and dry the impression thoroughly. Once the first impression has fully set take the second one with Doric Quicktime Flo-Light, Doric Flo-Light or Doric Light, applying it around the gingival margins and over the occlusal surfaces and the set putty. Seat the impression slowly so that any excess material can escape without building up internal stresses. When seated, hold the tray in position with a light, firm touch.
- When the Doric wash material has set, remove the tray from the mouth and proceed as normal.

Storage

Store in a cool dry place at a temperature below 23°C.

Health & Safety

In case of patients with a known allergy to polyethylene discontinue use and seek medical advice. Keep silicone materials away from eyes and clothing and in case of contact, wash thoroughly with soap and water. For further information request Health & Safety Data Sheet.

Lot numbers and expiry dates are shown on all containers.
See above for storage instructions.

This product is specifically formulated for use in dentistry.

Davis Schottlander & Davis Ltd.
Fifth Avenue, Letchworth, Herts SG6 2WD England
Tel:+(0)1462 480848 www.schottlander.com



W3570/0908