

Natura Denture Teeth Instructions for Use

INTRODUCTION

Natura are polymer artificial teeth for dental prostheses in accordance with EN ISO 22112:2017.

NOTE. The denture making process is undertaken by professionally qualified dentists, clinical dental technicians and dental technicians. These IFUs do not attempt to specify the work for which they are trained such as taking impressions of the mouth, identifying anatomical landmarks or the occlusal scheme to be employed since these processes are within their professional competencies and remit. Nor do they attempt to define the laboratory workflow since these are not functions of the teeth themselves. These instructions are applicable when the lost wax technique is employed for production of the denture base.

INDICATIONS

For use in the manufacture of full and partial dentures including implant-supported dentures.

CONTRA-INDICATIONS

Not recommended for the manufacture of fixed crowns and bridges.

SETTING UP THE TEETH

Prepare and mount models on an articulator in accordance with information provided by the clinician and standard laboratory practice. Unless already chosen by the clinician, choose the moulds of teeth to be used from the Natura mould chart or chairside assortment and the shade from the Natura shade guide. Once the upper anterior mould has been chosen, the mould chart makes suggestions as to which sizes of lower anterior and posterior moulds best fit that of the chosen upper anterior mould. This will depend on the clinical situation and personal preference. Lay a wax base onto the model, trim the teeth where necessary to fit over the ridge and mount them in accordance with the clinical information provided. Check articulation on an articulator.

FINAL WAXING UP

Final waxing up should be carried out paying particular attention to the anatomical details that will be highlighted in the final denture. Time spent at this stage will reduce the time needed in finishing.

WAX REMOVAL

When investing the denture in plaster, immerse the flask in boiling water for five minutes. Remove and open the flask, lift out wax and flush thoroughly with boiling water to which detergent solution has been added, followed by clean boiling water. Allow to cool.

APPLICATION OF ISOLATING SOLUTION AND INVESTING

Use a clean brush to apply one coat of Schottlander Isolating Solution to all plaster surfaces. Absorb excess from around teeth while it remains fluid. Make quite sure there is none on the exposed surfaces of the teeth. Allow to dry before investing which should be done in the normal manner.

FINISHING

The care taken with waxing and flasking will leave very little finishing to be done. Where plaster has been used place in an ultrasonic bath to remove all plaster from between the teeth. *CARE.* Plaster removal solutions which contain sodium hydroxide (caustic soda) may cause discolouration of both enigmalfite and other brands of composite denture teeth. Their use should therefore be avoided. Schottlander Plaster Removal Solution is a suitable product for use in this application. Trim away any flash and smooth and pumice the borders of the denture and palatal area (if desired). Limit polishing of the teeth to a light buffing. In these areas avoid the use of pumice which will tend to remove the outer surface detail of the teeth and any stippling of the denture base.

STORAGE

Store Natura teeth in a cool dry place and avoid prolonged exposure to sunlight. Keep boxes or cabinets closed when not in use.



HEALTH AND SAFETY

P261 avoid breathing dust/fumes/ gas/mist/vapours/spray.



LOT NUMBERS

The Lot number and use by dates are shown on all tooth cards.

This product is specifically formulated for use in dentistry.

Natura is an internationally registered trademark of Davis Schottlander & Davis Limited.

Where this product is shown as having been certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, this is exclusively for the indication(s) shown in the above Instructions for Use. Other non-medical uses ascribed to this device are not within the scope of CE certification, and users should be aware that product performance and/or safety has not been evaluated by SGS for those purposes.

 **Schottlander (IRL) Ltd.**
Suite 11, Parklands Office Park, Southern
Cross Road, Bray, Co Wicklow, Eire.

 **UK
CA**
0120

 **MD**

 **CE**
1639