

R-Motion

INFORMATION CARD



Reciprocation Redefined



R-Motion® Glider
Reciprocating glide path instrument




R-Motion® shaping file
R-Motion 30 (30/.04)
or R-Motion 25 (25/.06)



- Create straight-line coronal and radicular access
- Use an ISO 010 hand file to explore the canal
- Determine the working length (WL)
- Progress with **R-Motion Glider** with **gentle 2-3mm strokes** (max. 3), applying very light apical pressure and allowing file to passively progress along the canal until reaching the WL
- Irrigate & recapitulate
- Use the appropriate **R-Motion shaping file** with **gentle 2-3 mm strokes** (max. 3), applying very light apical pressure and allowing file to passively progress along the canal until reaching the WL
- Irrigate & recapitulate
- Once the canal has been shaped, proceed with the cleaning of the canal
- Obturate by selecting the appropriate GP point and sealer

**GLIDE
PATH FILE**

SHAPING FILE SELECTION*



This product is specifically manufactured for use in dentistry.

Select the final shaping file depending on the canal size and your own therapeutic approach.

*Select the **R-Motion** single-file shaping instrument which is best suited to the anatomy of the canal being treated. If the radiographs do not allow you to clearly identify the anatomy of the canal, it can be gauged using hand files.