

Matchmaker Zirconia Framework Finishing Kit

Maximum Operating Speeds

LBZ-2 **Coarse**  3' - 18000 rpm

LBZ-3 **Coarse**  3' - 6000 rpm

LBZ-4 **Medium**  3' - 10000 rpm

LBZ-5 **Medium**  3' - 8000 rpm

LBZ-6 **Super Fine**  3' - 10000 rpm

LBD-3 **Medium**  - 300000 rpm

LBD-27 **Medium**  - 300000 rpm

LDD -7  3' - 5000 rpm