

Scout & Glide Path RaCe

Rotary Nickel Titanium Files for special situations

START POINT

Scout RaCe - use in canals with severe curvatures where it is difficult to get a hand file larger than #6 or 8 to working length.

Severe Curvatures

Glide Path RaCe - use in highly sclerosed canals for coronal/middle third enlargement.

Sclerosed Canals

Scout RaCe Kit
Contains 6 instruments: 2 each of .02 10, 15 & 20
Code: 849-5

Scout RaCe instruments in individual packs of 5
Taper .02

Taper	ISO	25mm
.02	10	837-50
.02	15	837-51
.02	20	837-52

Glide Path RaCe instruments in individual packs of 5
Tapers .02, .04 & .06

Taper	ISO	25mm
.02	10	837-50
.04	10	836-50
.06	10	835-50

Davis Schottlander & Davis Ltd.
Fifth Avenue, Letchworth Garden City, Hertfordshire SG6 2WD England
Telephone +44(0)1462 480848 www.schottlander.com

Freephone 0800 97 000 79

This product is specifically manufactured for use in dentistry. 0120

Products and training for better dentistry

schottlander

WG3670/0711