

Dr Christopher Orr Posterior Restoration & Equilibration Kit

Product
Code 188-P

INFORMATION CARD

**Direct
Restoration
Preparation**

Note

Use

Quick Order Code



Round End
Cylinder

Caries
access

M296



Tungsten CX
Cross Cut

Remove old
amalgam

2158

**Inlay/Onlay
Preparation**

Note

Use

Quick Order Code



Tungsten CX
Cross Cut

Remove old
amalgam

2158



Round End
Taper

Internal
preparation

M313



Long Fine
Flame

Clear
contacts

F395



Rugby
Ball

Cuspal
reduction

M052



Rugby
Ball

Refine
preparation

F052

This product is
specifically
manufactured for
use in dentistry.

**Crown
Preparation**

Note

Use

Quick Order Code



Fine
Flame

Break
contacts

M385



Round End
Cylinder

Axial
preparation

M410



Round End
Cylinder

Axial
preparation

M240



Round End
Cylinder

Axial
preparation

M244



Occlusal

Occlusal surface
preparation

M041



Rugby Ball

Occlusal surface
preparation

M052



Rugby Ball

Refine
preparation

F052



End Cutting

Refine
margins

M360

**Other
Applications**

Note

Use

Quick Order Code



Equilibration
Burr

Adjust
ceramic/
occlusal
equilibration

M262



Rugby Ball

Adjust
ceramic/
provisionals

F052



Rugby Ball

Fine
adjustment
of ceramic

XF052